

SURGICAL SUPPLIES

SALE ON HADECO ULTRASONIC DOPPLERS



BIDOP ES100VS



HADECO BIDOP BI-DIRECTION WAVEFORM LCD VASCULAR DOPPLER WITH 8MHZ PROBE

The Bidop is a uniquely designed bi-directional pocket Doppler with LCD display, USB for optional Computer interface and optional PPG ability. It is supplied with an 8MHz Probe.

Features & Specification:

- LCD display 128 x 64 dots, STN LCD
- Bi-directional wave form (normal & slow mode) display
- Waveform Memory - 30 Waveforms
- Mode settings: Memory, Waveform, Direction, Time Scale & others
- Numerical Data Display: systolic, diastolic & mean Velocities.
- Speaker Output: 200mW
- Power - DC 9 Volts, Alkaline
- Battery Life: 2.5 hours
- Automatic shut off
- Size: 75 (W) x 140 (D) x 25 (H) mm
- Weight: 350 grams (including battery & probe)

NOW \$1290.00

EQ0037

HADECO MINIDOP VASCULAR DOPPLER WITH 8MHZ PROBE

This simple Doppler has great sounds and can be used with five interchangeable probes having different frequencies of 2, 4, 5, 8 and 10 MHz, to detect fetal heartbeats as well as arterial and venous blood flows in the extremities. Supplied with 8MHz Probe.

Features & Specification:

- DC9V, A9 volt alkaline, battery
- Weight : 270 grams.
- Multi freq :2,4,5,8,10 MHz
- Speaker Output 200mW.
- Dimension : 75 (w) x 140 (d) x25 (h) mm
- Accessories: Carrying case, Batteries, Earphones, 60ml Gel

MINIDOP ES100VX



NOW \$695.00

EQ0032

PODIATRY

All prices quoted without GST