

SURGICAL SUPPLIES

CLEARANAIL PRO FUNGAL DEVICE



Clearanail Pro is an innovative method of treating fungal nails. Hitting the UK with great success this computer controlled device painlessly drills microscopic holes into the infected nail. The nail is then sprayed with a powerful antifungal on a daily basis for 2-3 months. Because of the microscopic holes the active antifungal ingredient is able to permeate the nail and underlying nail bed where the infection is located. Advantages with this next generation model 50% faster drilling, up to 90 minutes per treatment.

Includes 5 Single Use Microcutter sets (SUM)
5 Resolve Solutions 25ml


\$6500.00



OVER 500 TREATMENTS ALREADY
CARRIED OUT IN NEW ZEALAND

PODIATRY

DIAFACTORY TINEA UNGUIUM DETECTION KIT



- *QUICK AND EASY TO USE
- *FAST RESULTS WITH 99% ACCURACY WITHIN 5 MINUTES
- *DETECTS IF ANY FUNGAL INFECTION PRESENT
- *AFFORDABLE TO PATIENTS AND CARRIED OUT IN CLINIC

\$250

BOX OF 10

All prices quoted without GST