

SURGICAL SUPPLIES

PRICES DROPPED ON NSK DRILLS



THE PODIASPRAY 40,000RPM SPRAY DRILL WITH LED LIGHT

EQ0124

- Speed adjustable from 5,000 to 40,000 rpm
- LED Lighting for illuminating the work area
- Separate on/off switch for the spray function
- Water tank with a capacity of 125 ml
- Dimensions: height 10 cm, depth 23.7 cm, width 24.8 cm. weight: 4.3 kg



HOT OFFERS

\$2580.00

Normal Price \$4200



MASTER VAC DRILL WITH PDV30 HANDPIECE

EQ0123

- Speed adjustment via the membrane keypad
- Clockwise / counter-clockwise rotation
- Suction power adjustable via membrane keypad
- Suction power display - Suction power up to max. 90 mbar
- Noise level 35 to 59dB
- Dimensions: height 11.5 cm, width 28.5 cm, depth 24 cm. weight 3.4 kg

\$1680.00

Normal Price \$2710

THE PODIASPRAY CAP 30,000RPM SPRAY DRILL

EQ0156

- Powerful compact brushless motor
- Speed adjustable from 3,000 to 30,000 rpm
- Large easy-to-read speed display
- Clockwise/ counter-clockwise rotation
- On/off-switching of the spray function
- Individual spray adjustment
- Water tank with a capacity of 250 ml
- Dimensions: height 18.4 cm, depth 19.5 cm, width 27.3 cm. weight: 3.8 kg



\$1680.00

Normal Price \$2700

PODIATRY